


GLUCOSE DEHYDROGENASE

(FAD-DEPENDENT)

REACTION:



PRODUCT DESCRIPTION

Catalog No.:	SDZ500420
Appearance:	Yellow amorphous powder
Source:	Microorganism
Enzyme Commission Number:	EC 1.1.99.10
CAS Number:	37250-84-3
Storage temperature:	-20°C
Specific activity:	≥ 800U/mg protein
Unit definition:	One unit will oxidate one micromole of glucose per minute at pH 6.5 at 37°C.

PROPERTIES

Molecular weight:	97kD (SDS-PAGE)	
Isoelectric point:	6.0	
Michaelis constant:	6.0 × 10 ⁻² M (α D-Glucose)	
Optimum pH:	7.0	{Fig. 1}
Optimum temperature:	55°C	{Fig. 3}
pH Stability:	5.0-9.0 (25°C, 20hr)	{Fig. 2}
Thermal stability:	< 40°C (pH7.0, 30min)	{Fig. 4}
Inhibitors:	Cu ²⁺ , Zn ²⁺	
Effect of various chemicals:		{Table 1}

Table 1.**Effect of Various Chemicals on GDH (FAD-dependent)**

[The enzyme dissolved in 50mM Tris-HCl buffer with 0.1% BSA, pH7.5 (10U/ml) was incubated with each chemical at 37°C for 1.5hr.]

Chemical	Concn. (mM)	Residual activity
None	-	100%
CaCl ₂	2.0	99%
CoCl ₂	2.0	97%
CuSO ₄	2.0	77%
FeCl ₃	2.0	99%
MgSO ₄	2.0	100%
MnSO ₄	2.0	105%
NiCl ₂	2.0	83%
ZnSO ₄	2.0	55%

Chemical	Concn. (mM)	Residual activity
BME	2.0	95%
NEM	2.0	89%
EDTA	5.0	92%
NaN ₃	20.0	91%
Na-cholate	0.10%	99%
SDS	0.05%	94%
Triton X-100	0.10%	105%
Tween 20	0.10%	104%
Boric Acid-Borax	2.0	101%

