


ACETYLPIRIDINE-ADENINE DINUCLEOTIDE

PRODUCT DESCRIPTION

Cat. No.:	SDZ600080
Product:	3-acetylpyridine-adenine dinucleotide ; 3-acetylpyridine-adenine dinucleotide, reduced form; 3-acetylpyridine nad;
Appearance:	Light yellow powder
CAS Number:	86-08-8
Solubility:	Soluble in water
Storage temperature:	<20 °C
Formula:	$C_{22}H_{28}N_6O_{14}P_2$
Formular weight:	662.4
Purity:	$\geq 84\%$ (by HPLC)



APAD