


# **PHOSPHOGLUCOMUTASE**

α-D-glucose 1,6-phosphomutase

## **REACTION:**



## **PRODUCT DESCRIPTION**

Catalog No.:	SDZ500090
Appearance:	White power, lyophilized
Source:	Rabbit muscle
Enzyme Commission Number:	EC 5.4.2.2
CAS Number:	9001-81-4
Storage temperature:	-20℃
Specific activity:	≥150U/mg protein
Unit definition:	One unit will convert one micromole of α -D-glucose 1-phosphate to α -D-glucose 6-phosphate per minute at pH 7.4 at 30℃.

**PGM**